



DEFINING A PROCESS FOR STRESS REDUCTION IN THE KEEL TRAY INTERFACE IN UNICOMPARTMENTAL KNEE REPLACEMENT TIBIAL COMPONENTS

Laurence Marks, Stephen Mellon, David Murray

Oxford Orthopaedic Engineering Centre

Nuffield Department of Orthopaedics, Rheumatology and Musculoskeletal Sciences

University of Oxford



The Oxford Unicompartmental Knee replacement

- 3 Part device
 - Femoral component
 - Meniscal bearing
 - Tibial component
- Design improvement of the tibial section using finite element non-parametric shape optimisation.




Design Driver

- Minimum bone resection is advantageous
- This can be achieved through reducing the tray thickness



Tray Failure

- Structural performance and durability constrain geometric design freedom
- Hasn't ever happened to an Oxford Unicompartmental Knee, but has happened with other devices.






Initial experience of the Journey-Deuce bicompartamental knee prosthesis. Palumbo BT, Henderson ER, Edwards PK, Burris RB, Gutiérrez S, Raterman SJ. *J Arthroplasty* 2011; 26 (6) Suppl 1: 40-45.

ASTM testing




- Test specification

- Test vs real life vs simulation trichotomy
- ASTM test doesn't replicate real world usage
- Almost into codes based design here
- But its probably the best test case we can apply..

- Finite element models of the test can provide design insights




Full block support








Simulating the ASTM test

- CAD models: NX and Solid Edge
- FEA models: Abaqus CAE and solver.
- Geometry:
 - Implant geometry used without simplification
 - Bearing and Femoral component – used to apply load in realistic manner
 - Rollers – discrete rigid surfaces
 - Femoral device – reduced geometry, discrete rigid surface
- Contact:
 - Surface to surface contact
- Higher order tetrahedral elements in all deformable components.
- High levels of local mesh refinement



Simulating the ASTM test

- A range of sizes were analysed - stress results were always dominated by the stress in the keel tray intersection.
- We didn't have any fatigue allowable's for the material, or knock down factors for surface finish etc
- So in this study we used the results of a device know to pass the test as a benchmark
- But this is a workflow presentation on design improvement processes..




Simulating the ASTM test

- Challenges

- Solution timescales – typically 4-24hrs depending on model geometry, contact convergence etc.
- Instability – loading protocol developed to control

- BCs (5)
 - + fix_rollers
 - + initial_load
 - + restrain_bearing
 - + restrain_device
 - + restraint_for_forces



Appendix 1 – Load application protocol


As the number of contact interfaces and components means that the problem is very unstable a 3 stage loading approach has been adopted.

Step 1: Apply initial load using a small vertical displacement (0.05mm or 0.1mm). All rotations of load application point are fixed. 3 points on implant and 3 points on bearing fully restrained.

Step 2: Settle – restraints on implant and bearing removed. 4 steps run to allow new load system to equilibrate (the stresses at the restraint points on the device seem to take a whole step to die away.) Record reaction at supports at the end of this step.


Step 3: Load application – apply the load to the reference point on the cylinder. This is done using a unit load applied to the cylinder reference point and an amplitude which goes from the measured reaction at the end of step 2 to the peak load.

Step 4: Check reactions





The need for a rapid solution

- Design improvement means we need to investigate the designspace, but each point is computationally very expensive.
- Reduce model complexity to increase solution speed with minimal reduction in solution accuracy.
 - Replace contact regions with line constraints or load patch (size taken from contact pressure region on contact model).
 - Redraw simplified geometric representation – not simplify original.
- Solution times fall from many hours to several minutes.






The rapid solution?

- Less than 10% difference between complex contact model and simplified case
- Changing model geometries is rapid and straightforward. (And enables automation)





Looking at the design envelope

- We need a reduced stress raiser but we've got a restricted space envelope




Optimising the fillet radius

- Bigger fillet radius works, but exceeds the design restriction region
- Non-parametric shape optimisation offers a way forwards.




Some solve time metrics

- Multicore solution of contact model: min 4hrs
- Number of solution cycles for optimisation: approx. 60
- Number of optimisation cycles to achieve useful solution: typically 4
- Which gives us 960 hours solve time
- **40 days**
- It's a sequential process so cloud, cluster, parallel, GPU won't help you.
- Which effectively means this is unsolvable without a simplified model




Implementation of the process Abaqus /TOSCA NDORMS



Results of non-parametric optimisation


- Non-parametric shape gives significant stress reduction and significant reduction in stress gradient

2mm to 0.5mm fillet




Peak – 142MPa

0.5mm fillet






Peak – 250MPa

1mm fillet




Peak – 175MPa




Transferring the shape data

- To be useful we have to convert the non-parametric nodal positions into parametrically defined forms which can be scaled across the range of device sizes.
- Defining a CAD model form that didn't re-introduce stress raisers proved challenging




Lack of tangency in this region has reintroduced a stress raiser.



“Validate” simple model using complex interaction model


- Optimised design analysed both ways..



Tray thickness results


- New fillet geometry reduces stress across a range of tray thicknesses.
- Critically allowing a reduction of tray thickness without increase in stress level.

Max principal stress MPa



Process Map

- Mapped process reduces optimisation timescales and effort to industrially relevant levels.



Conclusion

- Non-parametric shape optimisation of the critical radius has unlocked regions of the designspace, improving component performance. The new fillet radius has allowed us to create a safe design with a thinner tray.
- This technology (and other optimisation and design space exploration techniques) can only be applied when simulation models are optimised and reduced in order to create meaningful, industrially relevant, analysis timescales.

